

# Copper - Gold Concentrate

## DIGESTION OF COPPER – GOLD CONCENTRATE USING COLDBLOCK™ DIGESTION TECHNOLOGY

*Author: Lorie-Anne Doig and Erick Helmeczi*

### Introduction

This application note will focus on the digestion of Copper – Gold concentrate using ColdBlock™ Digestion CB12L Technology.

### Method

A sample of CRM (Certified Reference Material) OREAS-991 was digested using the following two methods:

- 0.25g of sample was weighed and placed into a ColdBlock™ Digestion vessel.
- 12ml inverse Aqua Regia (1:3 HCL:HNO<sub>3</sub>) was added.
- Sample was digested at 65% power for 15 minutes.
- Chiller temperature was set to -5°C.
- Sample was allowed to cool.
- Sample was normalized to 50 mL with 20% HCl.

or

- 30g of sample was weighed and placed into a ColdBlock™ Digestion vessel.
- 160mL inverse Aqua Regia (1:3 HCL:HNO<sub>3</sub>) was added.
- Sample was digested at 65% power for 15 minutes.
- Chiller temperature was set to -5°C.
- Sample was allowed to cool.
- Sample was normalized to 200mL with 20% HCl.

### Instrument

ColdBlock™ Digestion CB12L Technology.

### General

This procedure is specific for the sample digested and may need modification for different samples to achieve the desired result.