

IARM-59A – Incoloy 825

SUMMARY

This application note is for the digestion of IARM 59A, Incoloy Alloy 825.

Instrument: ColdBlock CBM, ICP-MS, ICP-OES

Published: November 2023

Digestion Time: 15 Minutes

Acid Used: Aqua Regia, HF

Average ColdBlock Recovery vs. CRM:


- 101% Chromium
- 100% Iron
- 99% Molybdenum

METHODOLOGY

1. Chiller temperature was set to -5°C
2. 0.25g of IARM 59A was weighed and placed into a ColdBlock HF friendly test tube
3. 12mL of Aqua Regia + 1mL HF was added
4. Sample was digested at 90% power for 15 minutes
5. Sample was cooled and bulked to 50mL using 2% HNO_3 v/v

DISCUSSION

- Samples were digested triplicate
- After 15 minute digestion, the samples are green and clear
- As a safer alternative to HF, you can add solid NH_4HF_2 (Ammonium Bi-fluoride) or NH_4F (Ammonium Fluoride)
- IARM 59A is in the form of chips



IARM 59A after topping up to 50mL

IARM 59A- Inco 825 is a certified reference material (CRM) sourced from LGC ARMI. LGC ARMI, Analytical Reference Materials International, Manchester, New Hampshire, USA (October 1992)

IARM-59A – Incoloy 825

Results

IARM 59A								
Analyte:	Certified Value (wt. %)	+/-	ColdBlock 1 (wt. %)	ColdBlock 2 (wt. %)	ColdBlock 3 (wt. %)	Average	% RSD	Recovery
Al	0.05	0.0063	0.05	0.045	0.047	0.05	3.9	94%
Co	0.25	0.0133	0.248	0.259	0.247	0.25	2.2	101%
Cr	22.12	0.054	22.26	22.39	22.68	22.44	0.8	101%
Cu	1.71	0.0216	1.73	1.7	1.71	1.71	0.7	100%
Fe	30.82	0.0439	30.67	30.76	30.91	30.78	0.3	100%
Mn	0.33	0.01	0.32	0.31	0.31	0.31	1.5	95%
Mo	2.68	0.02	2.65	2.57	2.73	2.65	2.5	99%
Ni	40.91	N/A	40.66	41.1	41.1	40.96	0.5	100%
Nb	0.02	0.004	0.0208	0.0209	0.0207	0.02	0.4	104%
Si	0.1	0.0112	0.12	0.1	0.11	0.1	8.1	109%
Ti	0.83	0.0111	0.82	0.82	0.83	0.82	0.4	99%
V	0.03	0.0096	0.029	0.028	0.029	0.03	1.3	96%
W	0.13	0.013	0.13	0.13	0.13	0.13	0.3	101%
Zr	<0.01	N/A	0.00032	0.00038	0.00027	0.0003	13.9	N/A