

HIGH GRADE TUNGSTEN-MAGNETITE

DIGESTION OF CRM OREAS 701, A HIGH-GRADE TUNGSTEN-MAGNETITE ORE USING COLDBLOCK™ DIGESTION TECHNOLOGY

Author: Lorie-Anne Doig

Introduction

This application note will focus on the digestion of CRM OREAS 701 using ColdBlock™ Digestion Technology.

Method

Six samples of OREAS-701 were digested using the following method:

- Chiller temperature was set to -5°C.
- 0.25g sample was weighed and placed into a ColdBlock™ Digestion vessel.
- 6mL Phosphoric acid was added (H₃PO₄)
- Samples were digested at 85% power for 10 minutes.
- Samples were cooled and diluted to 50mL using 1% HNO₃.

Instrument

ColdBlock™ Digestion CB6L Technology.

General

This procedure is specific for the sample digested and may need modification for different samples to achieve the desired result.