

Gold Ore

DIGESTION OF A GOLD ORE USING COLDBLOCK™ DIGESTION TECHNOLOGY

Author: Lorie-Anne Doig

Introduction

This application note will focus on the digestion of a Gold Ore using ColdBlock™ Digestion CB12L Technology.

Method

A sample of Gold Ore was digested using the following method:

- 30g of sample was weighed and placed into a ColdBlock™ Digestion vessel.
- 140ml of inverse Aqua Regia (1:3 HCL:HNO₃) was added.
- Sample was digested at 100% power for 15min.
- The chiller was set to -5°C.
- Sample was allowed to cool.
- Sample was normalized to 200mL using 2% HCl + 1% L-cysteine

Instrument

ColdBlock™ Digestion CB12L Technology.

General

For samples rich in Sulphur add HNO₃ slowly first.

For samples high in silver, normalize to 200ml using 20% HCl + 1% L-cysteine.

This procedure is specific for the sample digested and may need modification for different samples to achieve the desired result.

