

## APPLICATION NOTE

# Estuarine Sediment, NIST 1646a

## SUMMARY

The application note summarizes the digestion of NIST 1646a, an Estuarine Sediment standard reference material using ColdBlock Digestion Pro Series Technology.

<b>Instrument:</b>	ColdBlock CBM, chiller, ICP-OES, ICP-MS (with CRC technology)
--------------------	---

<b>Published:</b>	March 2025
-------------------	------------

<b>Digestion Time:</b>	45 Minutes
------------------------	------------

<b>Acid Used:</b>	HNO <sub>3</sub> , HCl & H <sub>2</sub> O <sub>2</sub>
-------------------	--

<b>Average ColdBlock Recovery vs. CRM:</b>	<ul style="list-style-type: none"><li>■ 102% Barium</li><li>■ 107% Chromium</li><li>■ 102% Molybdenum</li></ul>
--	---

## METHODOLOGY

1. Chiller temperature was set to -5°C
2. 0.5g of each sample was weighed and placed into a ColdBlock Digestion vessel
3. 10mL of 1:1 HNO<sub>3</sub> + 5mL 1:1 HCl were added, and the samples were digested at 70% power for 5 minutes
4. 10mL HNO<sub>3</sub> was added and digested again at 60% power for 20 minutes (repeat this step if fuming continues)
5. 3mL of 30% H<sub>2</sub>O<sub>2</sub> was added slowly, and digested again at 60% for 10 minutes
6. 5mL HCl was added and digested at 60% power for 10 minutes
7. Samples were cooled and bulked to 50mL using ultrapure water



15 samples of NIST 1646a after bulk up to 50mL

## DISCUSSION

- This method yielded data comparable to EPA 3050B standards in just 45 minutes
- The addition of HCl in the final step of the digestion helps keep certain elements such as antimony and silver solubilized
- This method is not appropriate for standard ICP-MS instruments that are not equipped with Collision/Reaction Cell technology

## Results

Estuarine Sediment, NIST 1646a					
Method:	0.5g - Add 10mL 1:1 HNO <sub>3</sub> + 5mL 1:1 HCl, digest at 70% for 5 minutes, then add 10mL HNO <sub>3</sub> and digest at 60% for 20 minutes (repeat this step if fuming continues) Add 3mL of 30% H <sub>2</sub> O <sub>2</sub> , slowly and digest again @ 60% for 10 minutes, then add 5mL HCl and digest again @ 60% for 10 minutes. Cool and bulk to 50mL with ultrapure water				
Element	*Reference Value (mg/kg) n=5	ColdBlock Average (mg/kg) n=16	ColdBlock Stdev	ColdBlock % RSD	ColdBlock % Recovery vs SW-846 3050B
Aluminum	7075	7182	206	2.9%	102%
Antimony	0.72	0.81	0.04	5.1%	112%
Arsenic	4.69	6.12	0.28	4.5%	130%
Barium	21.2	21.7	1.03	4.7%	102%
Cadmium	0.127	0.142	0.00	2%	112%
Calcium	3578	3711	162	4.4%	104%
Chromium	24.4	26.2	2.15	8.2%	107%
Cobalt	3.2	3.4	0.13	3.8%	108%
Copper	9.98	9.96	0.13	1.3%	100%
Iron	15330	17111	849	5.0%	112%
Lead	5.26	4.99	0.3	5.1%	95%
Magnesium	2972	3347	35	1.1%	113%
Manganese	117.3	125.2	1.7	1.3%	107%
Mercury	0.03	0.05	0.0039	8.5%	142%
Molybdenum	1.81	1.84	0.15	8.4%	102%
Nickel	18.8	21.9	2.27	10.4%	116%
Potassium	2235	2273	96	4.2%	102%
Sodium	3819	3913	254	6.5%	102%
Vanadium	20.95	22.90	1.10	4.8%	109%
Zinc	36.72	35.92	2.06	5.7%	98%