# Geostats - GCr-06 - Chromite

#### **SUMMARY**

The application note summarizes the digestion of GCr-06, a Chromite pulp Certified Reference Material using ColdBlock™ Digestion Pro Series Technology.

Instrument:	ColdBlock CBM sample digester, chiller, Quartz vessels, ICP-OES				
Published:	February 2023				
Digestion Time:	20 Minutes				
Acid Used:	$H_3PO_4 \& H_2SO_4$				
Average ColdBlock Recovery vs. CRM:	<ul><li>99% Chromium</li><li>99% Iron</li><li>100% Magnesium</li></ul>				

#### **METHODOLOGY**

- 1. Chiller temperature was set to -5°C
- 0.25g of sample was weighed and placed into a ColdBlock™ Digestion vessel
- 3.  $8 \text{mL of H}_3 \text{PO}_4 \& 4 \text{mL H}_2 \text{SO}_4 \text{ was added}$
- 4. Sample was digested at 80% power for 20 minutes
- 5. Samples were cooled and bulked to 100 mL using 2% HNO<sub>3</sub>

#### **DISCUSSION**

- Phosphoric acid can etch borosilicate glass so ColdBlock's Quartz test tubes were used
- After 20 minutes, samples are green in color and no solid material remains
- Samples were analyzed by ICP-OES



Figure 1 - GCr-06 after bulk-up to 100mL

Prior to homogenization and testing, this material was sourced from Pulp Chromite  $\,$ 

geostats.com.au

GCr-06; Pulp Chromote; Geostats Pty Ltd, Mining Industry Consultants; O'Connor, Western Australia (December, 2010)

## Geostats - GCr-06 - Chromite

### **Results**

Geostats - GCr-06 - Chromite Pulp												
Method:	0.25g	8mL H <sub>3</sub> PO <sub>4</sub> - Digested at 80% power for 20 minutes										
Element	Certified XRF Value (%)	95% Confidence Limits								%		
		Low	High	Sample A	Sample B	Sample C	Average (ppm)	Stdev	% RSD	Recovery vs certified XRF		
Al <sub>2</sub> O <sub>3</sub>	8.72	8.68	8.76	8.68	8.69	8.68	8.68	0.005	0.1	100%		
Cr <sub>2</sub> O <sub>3</sub>	47.92	47.72	48.12	47.61	47.55	47.64	47.60	0.037	0.1	100%		
Fe <sub>2</sub> O <sub>3</sub>	28.8	28.7	28.9	28.74	27.79	29.02	28.52	0.526	1.8	99%		
MgO	9.01	8.96	9.06	8.92	8.99	9.01	8.97	0.039	0.4	100%		
MnO	0.716	0.709	0.723	0.734	0.745	0.729	0.736	0.007	0.9	103%		
TiO <sub>2</sub>	0.343	0.337	0.349	0.329	0.339	0.341	0.336	0.005	1.6	98%		