

APPLICATION NOTE

Copperbelt, Sulphide ZM - AMIS 0809

SUMMARY

The application note summarizes the digestion of AMIS 0809, a Copperbelt sulphide ore from Zambia, using ColdBlock™ Digestion Pro Series Technology.

Instrument: ColdBlock CBM Pro Series, chiller, ICP-OES

Published: June 2025

Digestion Time: 35 Minutes

Acid Used: HNO_3 , HCl , HF & H_3BO_3

Average ColdBlock Recovery vs. CRM:


- 105% cobalt
- 103% copper
- 102% nickel

METHODOLOGY

1. Set the chiller temperature to -5°C .
2. Weigh 0.25 g of each sample and transfer into a ColdBlock™ Digestion vessel
3. Add 20 mL of reverse aqua regia (add HNO_3 first, followed by HCl) and 3 mL of HF to each vessel.
4. Digest samples at 80% power for 20 minutes.
5. Add 20 mL of 4% boric acid (w/v) to each sample.
6. Digest samples again at 80% power for 10 minutes.
7. Add 5 mL of HCl .
8. Digest at 100% power for 5 minutes.
9. Cool the samples and adjust the volume to 50 mL with DI water

DISCUSSION

- The addition of boric acid aids in the re-solubilization of insoluble fluorides and neutralizes residual hydrofluoric acid (HF). This step is crucial for sample dissolution and minimizing analyte loss during digestion.
- Following the digestion process, the samples were visibly clear, indicating that the majority sample matrix had dissolved. A trace amount of sample material settled on the bottom of the tube overnight.



AMIS 0809 was made using sulphide ore sourced from the Kansanshi Mine, located in the Northwestern Province of Zambia. African Mineral Standards (AMIS). (2022). Certificate of Analysis: AMIS0809 - Copperbelt, Sulphide Zambia. Retrieved from www.amis.co.za/

Results

Copperbelt, Sulphide ZM - AMIS 0809

Method:	0.25g - Add 20mL reverse aqua regia + 3mL HF and digest at 80% power for 20 minutes. Add 20mL of 4% boric acid _{w/v} and digest again at 80% power for 10 minutes. Add 5mL HCl, and digest again for 5 minutes. Let cool, and adjust the volume to 50mL with DI water.				
Element	AMIS Certified Values	ColdBlock Average	ColdBlock Stdev +/-	ColdBlock RSD	ColdBlock Recovery vs Certified Value
Ag (ppm)	0.6	<1	N/A	N/A	N/A
Al (%)	1.04	1.00	0.02	1.8%	96%
As (ppm)	32	40	1.51	3.5%	114%
Ba (ppm)	21	20.6	0.44	2.1%	98%
Be (ppm)	0.2	<0.5	N/A	N/A	N/A
Bi (ppm)	1	<5	N/A	N/A	N/A
Ca (%)	18.15	17.97	0.21	1.2%	99%
Cd (ppm)	0.2	<0.5	N/A	N/A	N/A
Co (ppm)	366	387	1.52	0.4%	105%
Cr (ppm)	99	98	2.41	2.5%	99%
Cu (%)	2.97	3.05	0.02	0.7%	103%
Fe (%)	5.69	5.94	0.08	1.4%	104%
Mg (%)	1.22	1.23	0.01	1.1%	101%
Mn (ppm)	1065	1065	8.61	0.8%	100%
Mo (ppm)	100	108	3.69	3.4%	107%
Na (ppm)	8796	8686	183	2.1%	99%
Ni (%)	310	317	5.23	1.6%	102%
P (ppm)	291	286	11.9	4.2%	98%
Pb (ppm)	29	30	3.23	10.1%	101%
S (%)	5.29	5.52	0.07	1.3%	104%
Sb (ppm)	2	<5	N/A	N/A	N/A
Sr (ppm)	174	167	0.84	0.5%	96%
Zn (ppm)	75	71	0.67	0.9%	95%