

# Lead Concentrate, CPb-3

## SUMMARY

The application note summarizes the digestion of CPb-3, a lead concentrate certified reference material using ColdBlock Digestion Pro Series Technology.

**Instrument:** ColdBlock CBM Pro-Series, chiller, HF liners, ICP-OES

**Published:** May 2025

**Digestion Time:** 30 Minutes

**Acid Used:**  $\text{HNO}_3$ , HCl, HF &  $\text{H}_3\text{BO}_3$

**Average ColdBlock Recovery vs. CRM:**

- 100% antimony
- 102% lead
- 100% zinc

## METHODOLOGY

1. Set the chiller temperature to  $-5^\circ\text{C}$ .
2. Weigh 0.25 g of each sample and transfer into a ColdBlock Digestion vessel with liner
3. Add 20 mL of reverse aqua regia (add  $\text{HNO}_3$  first, followed by HCl) and 3 mL of HF to each vessel.
4. Digest samples at 80% power for 20 minutes.
5. Add 20 mL of 4% boric acid ( $w/v$ ) to each sample.
6. Digest samples again at 80% power for 10 minutes.
7. Cool the samples and adjust the volume to 50 mL with DI water




Figure 1- CPb-3 after bulking-up to 50mL

## DISCUSSION

- The addition of boric acid aids in the re-solubilization of insoluble fluorides and neutralizes residual hydrofluoric acid (HF). This step is crucial for sample dissolution and minimizing analyte loss during digestion.
- Following the digestion process, the samples were visibly clear - indicating that the sample matrix had dissolved. A trace amount of material settled on the bottom of the tube - see figure 1

## Results

CPb-3, Lead Concentrate					
Method:	0.25g - Add 20mL reverse aqua regia + 3mL HF and digest at 80% power for 20 minutes. Add 20mL of 4% boric acid <sub>w/v</sub> and digest again at 80% power for 10 minutes. Let cool, and adjust the volume to 50mL with DI water				
Element	CPb-3 Certified Values	ColdBlock Average n=5	ColdBlock Stdev +/-	ColdBlock RSD	ColdBlock Recovery vs Certified Value
Al (%)	0.203	0.200	0.005	2.5%	99%
Ca (%)	0.059	0.053	0.003	5.7%	90%
Cd (%)	0.0652	0.0653	0.002	3.1%	100%
Co (µg/g)	13.6	13.7	1.4	10.2%	101%
Cr (%)	0.0102	0.0104	0.0002	1.9%	102%
Cu (%)	0.240	0.245	0.004	1.6%	102%
Fe (%)	8.45	8.51	0.08	0.9%	101%
Mg (%)	0.1062	0.1	0.0024	2.3%	100%
Pb (%)	58.53	59.02	0.39	0.7%	101%
S (%)	17.03	17.93	0.21	1.2%	105%
Sb (%)	0.5800	0.5778	0.008	1.4%	100%
Zn (%)	5.96	5.98	0.04	0.7%	100%

The material for CPb-3 was donated by a North American refinery.

CCRMP CanmetMINING NRCan (2020). Certificate of Analysis for CPb-3, Lead Concentrate <https://natural-resources.canada.ca/>